## For CEREC, STML (Super Translucent Multi Layered) & ST (Super Translucent)

### [ENGLISH] INSTRUCTIONS FOR USE


#### I. Introduction

KATANA Zirconia Block is a pre-sintered zirconia block, which is designed for CEREC milling systems. (Please refer to your milling system's technical instructions for correct machine operation.) KATANA Zirconia block STML has 2 sizes: (a: 15.3x:19.2x:20.2 mm and a:17.8x:19.2x:20.2 mm) and has 14 shade variations: (A1, A2, A3, A3.5, A4, B1, B2, B3, C1, C2, C3, D3, D2, D3, NW). KATANA Zirconia block ST has 1 size: (a:15.3x:19.2x:20.2 mm) and has 1 shade variation: (CL). KATANA Zirconia Block STML consists of 4 graded shade layers: whered.KATANA Zirconia Block ST is mono color.

#### II. Intended Use

KATANA Zirconia Block is used for the fabrication of the ceramic restorations (FCZ crowns, inlays, onlays and veneers).

#### III. Image of product [For CEREC]

The STML is a color to vary from the enamel layer to the body layer through two gradation layers in the x-axis direction.

#### IV. Composition

ZrO<sub>2</sub>, Y<sub>2</sub>O<sub>3</sub> etc.

#### V. Type and Class (ISO6872:2015)

#### VI. Physical Properties

Coefficient of Thermal Expansion (25-500°C (77-932°F)): 9.8 (±0.2)×10<sup>-6</sup>/K

#### VII. Directions for Use

- Take the block from the packaging and confirm that the block does not have a crack or other damage.
- Place the block into the milling machine; then begin the milling process following the milling systems technical instructions.
- After milling, remove the restorations from the block with a diamond bur etc.

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

- The margins should be prepared with a deep chamfer and rounded shoulders, with cutting edges and corners rounded to eliminate sharp preparation corners. The angle of the axial surface should be within the angle of 5 to 15 degrees.
- When preparing teeth, avoid the following: deep shoulders, J-margins, knife edges, serrated margins, non-tapered abutments, undercuts, guide grooves, the formation of retentive holes, and sharp corners.
- Keep the following thickness of this product for fabricating prosthetics:

Location & indication	Wall thickness
Veneer	0.4 mm or more
Anterior crown	0.8 mm or more
Posterior crown, Inlay and onlay	1.0 mm or more

- Choose a shade color that is brighter than the intended color for a thick restoration, as it may look duller depending on the thickness of the restorations.
- Do not use the other sintering furnace except CEREC SpeedFire. When sintering this product for the first time, confirm with the supplier of CEREC SpeedFire that it is compatible with this product.
- Check the coefficient of thermal expansion of the porcelain glaze in the manufacturer's technical instructions to confirm compatibility.

**Caution in conjunction with usage:**

- In STML, the carved seal side of "inLab" is the body layer; the opposite side is the enamel layer. Depending on the setting of the scanner or milling machine, the enamel layer and the body layer may be machined in the opposite direction. Confirm that the block is compatible before processing.
- Dry milling is recommended. If wet milling is used, the translucency of the block may be reduced.
- Do not use the block if there is a crack noticed after removing it from the package.
- If you find a crack in the restorations, do not use.

**Warning:**

- If the patient is hypersensitive to zirconia or any other components, this product must not be used.
- Only use this product to produce restorations for anterior crowns, posterior crowns, inlays/ onlays and veneers.

**Warning:**

- The patient or the dental professional demonstrates a hypersensitivity reaction, such as rash, dermatitis etc., discontinue use of the product and seek medical attention immediately.
- This product should not be used when malocclusion, clenching or bruxism conditions are present.
- When milling "C" cutting, grinding and polishing the restorations, use an approved dust mask and vacuum with air filter to protect your lungs from inhaling dust.
- When milling the block or cutting, grinding and polishing the restorations, use eye protection to prevent the dust from getting into your eyes. If the dust gets into your eyes, immediately rinse with copious amounts of water and consult a physician.
- Do not use for any purposes except for dental restoration. This product is for dental application only.
- Do not touch the items heated by the furnace with your bare hands.
- This product is particularly translucent; therefore, one should carefully consider its appropriate use for restorations that require a translucent surface, such as a metal abutment, or a very dark or discolored tooth preparation.

#### VIII. Remarks on Handling

**Contraindications:**

- If the patient is hypersensitive to zirconia or any other components, this product must not be used.
- Only use this product to produce restorations for anterior crowns, posterior crowns, inlays/ onlays and veneers.

**Warning:**

- If the patient or the dental professional demonstrates a hypersensitivity reaction, such as rash, dermatitis etc., discontinue use of the product and seek medical attention immediately.
- This product should not be used when malocclusion, clenching or bruxism conditions are present.
- When milling "C" cutting, grinding and polishing the restorations, use an approved dust mask and vacuum with air filter to protect your lungs from inhaling dust.
- When milling the block or cutting, grinding and polishing the restorations, use eye protection to prevent the dust from getting into your eyes. If the dust gets into your eyes, immediately rinse with copious amounts of water and consult a physician.
- Do not use for any purposes except for dental restoration. This product is for dental application only.
- Do not touch the items heated by the furnace with your bare hands.
- This product is particularly translucent; therefore, one should carefully consider its appropriate use for restorations that require a translucent surface, such as a metal abutment, or a very dark or discolored tooth preparation.

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